

Ex-Vivo Cell Therapy Platform for Immune Reprogramming in Autoimmunity

O'Rourke, Brian., PhD1, Nguyen, Sunny Sang., BS1, Igo, Peter., BS1, Tilles, Arno W., MD, MSBE1, Parekkadan, Biju., PhD1-23, Barcia, Rita., PhD1

Sentien Blotechnologies, Inc. 99 Hayden Ave, Suite 140, Lexington, MA 02421 ² Department of Blomedical Engineering, Rutgers University, Piscataway, New Jersey 08854 ³Center for Surgey, Imovation, and Bloengineering, Massachusetts General Hospital, Boston, MA 02114 ³Harvard Stem Cell Institute, Cambridge, MA 02138

The Sentien Approach: Bringing Blood to the MSCs

Bioactive molecules secreted by MSCs are the primary source of activity of these therapeutically promising cells. We have engineered a system to maximize delivery of exposed to MSCs in SBI-101 therapy from MSCs to circumvent the half-life issues that have hindered MSC transplantation. This system overcomes the dosing constraints of IV infusion and potentiates a broad range of biological responses unparalleled in single molecule therapeutics.

Product Concept: Sentien is developing a combination product (SBI-101) that integrates allogeneic MSCs within an extracorporeal, blood contacting device to 2. fundamentally change the administration route. Instead of bringing MSCs to the blood, our product brings blood to the MSCs.

2. MSCs secrete therapeutic molecules creating concentrated, sustained microenvironment in SBI-101

Blood is continuously

Background: SBI-101-01 trial

Sentien currently has an ongoing Phase I/II clinical trial investigating SBI-101 in Acute Kidney Injury (AKI): Double blind, randomized, controlled, study at 2 doses to establish safety and pharmacologic POC (NCT03015623). Mare etal, Ridary 10 Rep. 2016





Mini-SBI-101 R&D Setup with Healthy PBMCs



Sentien has developed a miniaturized bioreactor system (mini-SBI-101) that enables detailed study of MSC effects on blood cells in the R&D setting, MSCs are seeded onto microreactors which are then perused using a pump. This miniaturized reactor system enables detailed study of MSC effects on blood cells.





MSCs in bioreactor inhibit lymphocyte proliferation. CFSE labelled PBMCs were activated with PHA and IL-2 and perfused through the MSC reactor for 5 days. Lymphocyte doublings were significantly reduced in the presence of MSCs

T cell activation surface markers CD35 are reduced in the presence of MSCs in bioreactor.

R&D of Mini-SBI-101 with PBMCs from Autoimmune Donors



PBMCs from donors with either SLE or MS are responsive to MSCs within the Bioreactor. PBMCs were activated with PHA and IL-2 and perfused through the MSC reactor for 5 days. PBMCS were then collected to assess immunophenotypic trades. Decases in T-cell markers CD4 and CD8 were noted, as well as increases in B-cell marker CD-19. Further, changes in the cytokine milleu were assessed by ELISA comparing between cellular and accellular treatment groups. Decreases in priorhfammatory cytokine were noted following exposure to MSCs.

Preliminary data supports broadening clinical application into autoimmunity conditions including SLE and MS

Summary

> TNFa

¥ L-23

A treatment option which could balance multiple dysregulated cytokine levels would provide a significant advancement in the treatment of autoimmune diseases. Here we have shown the R&D scaled mini-SBI-101 simultaneously modulated multiple targets associated with autoimmune disease.

For more information regarding how SBI-101 may provide a therapeutic option for the regulation of cytokine storm associated with COVID-19, please see:

Poster #836 Clinical Evidence for Immune Reprogramming with Extracorporeal Mesenchymal Stromal Cell Therapy



SBI-101

3. SBI-101 modulates immune cells leading to a systemic anti-

inflammatory response

Ex-vivo MSC therapy using SBI-101 technology has promise for many clinical applications requiring systemic immunotherapy for tissue repair and regeneration.

· MSCs